

**CONVERSION TABLE FOR LOGMAR TO SNELLEN'S
EQUIVALENT**

LOGMAR	SNELLEN EQUIVALENT
1.0	6/60
0.9	6/48
0.8	6/38
0.7	6/30
0.6	6/24
0.5	6/19
0.4	6/15
0.3	6/12
0.2	6/9.5
0.1	6/7.5
0.0	6/6
-0.1	6/5
-0.2	6/4
-0.3	6/3

Nb Each line on the chart represents a change of 0.1 log unit in the acuity level with each letter having a value of 0.02 log unit.
Therefore a patient reading correctly all letters on a specific line will score the full 0.1 log unit. For every extra letter on subsequent lines the patient reads correctly the patient will score an extra -0.02